

Effective Date: 12/30/25

Role(s): Clinical/Provider

Tacrolimus method change

NorDx will be switching methods from an Abbott Architect® assay to a Roche Cobas® E Assay for Tacrolimus testing.

This change in method may cause a shift in results. Due to this shift, results will not trend with previous patient history after the switch is made.

With this change, there will be a wider reporting range:

- Old range: 2 – 60 ng/mL
- New range: 2 – 90 ng/mL

No results will be reported outside the new range.

The reference ranges, potentially toxic, and critical values, will stay the same:

- ≤3 months post renal transplant: 12.0-15.0 ng/mL
- >3 months post renal transplant: 8.0-12.0 ng/mL
- Potentially toxic: 15.1-20.0 ng/mL
- Toxic: Critical Value >20.0 ng/mL

Please refer to the NorDx test catalog at <https://nordx.testcatalog.org/> for any processing or specimen collection questions. Please reach out to NorDx Customer Solutions if you have any questions – 207-396-7830.