

**Effective Date: 12/30/25**

**Role(s): Clinical/Provider**

## Homocysteine method change

NorDx will be switching methods from an Abbott Architect® assay to a Roche Cobas® C Assay for Homocysteine testing.

This change in method will likely result in an observed increase in homocysteine over previous testing. Due to this shift, results will not trend with previous patient history after the switch is made.

With this change, the reporting range will be less broad:

- Old range: 1 – 500 µmol/L
- New range: 3 – 250 µmol/L

No results will be reported outside the new range.

A new reference range will also go into effect.

- Old range: <14 µmol/L
- New range: <15 µmol/L

Please refer to the NorDx test catalog at <https://nordx.testcatalog.org/> for any processing or specimen collection questions. Please reach out to NorDx Customer Solutions if you have any questions – 207-396-7830