

THROMBOTIC RISK PROFILE COMPONENTS								
TEST NAME	Label Text	SPECIAL PT. INSTRUCTIONS	COLLECTION TUBE	PROCESSING INSTRUCTION	STORAGE/ TRANSPORTATION TEMP	TESTING LOCATION	TESTING DEPARTMENT	DEPARTMENT TESTING SCHEDULE
FACTOR V LEIDEN	FVDNA	NONE	1 LAVENDER (EDTA)	NONE - WHOLE BLOOD ASSAYS	ROOM TEMP OR REFRIGERATE	NorDx - SCARBOROUGH	MOLECULAR PATHOLOGY	WEEKLY (Wednesday)
PROTHROMBIN 20210 MUTATION	PTDNA							
HOMOCYSTEINE, PLASMA	HOMCY	FASTING preferred	1 LAVENDER (EDTA)	SPIN, AND SEPARATE OFF PLASMA (1.0 mL)	REFRIGERATE	NorDx - SCARBOROUGH	CLINICAL LAB	DAILY (M-F)
LUPUS ANTICOAGULANT HEMATOCRIT (Part of LUPUS test)	LUPUS (LAHCT)	NONE	1 LAVENDER (EDTA)	NONE - WHOLE BLOOD ASSAYS	REFRIGERATE	NorDx - SCARBOROUGH	CLINICAL LAB	DAILY (M-F)
ANTITHROMBIN III	ATIII	NONE	3 or 4 BLUE TOP (CITRATE)	DOUBLE SPIN, AND SEPARATE OFF PLASMA (5-7 ALIQUOTS CONTAINING 1.0mL)	FROZEN IMMEDIATELY	NorDx - BRAMHALL	SPECIAL COAG LAB	VARIES - TOTAL TAT IS ~21 DAYS
FIBRINOGEN	FIBRI	NONE						
LUPUS ANTICOAGULANT SCREEN (Part of LUPUS test)	LUPUS	NONE						
THROMBIN TIME	TCT	NONE						
INR	INRT	AFFECTED BY COUMADIN THERAPY						
PARTIAL THROMBOPLASTIN LTIME	PTT	NONE						
PROTEIN C ACTIVITY	PCACT	AFFECTED BY COUMADIN THERAPY						
PROTEIN S ACTIVITY	PSACT	AFFECTED BY COUMADIN THERAPY	Mayo Laboratories	Referral Testing	Weekdays Plus transit time			
BETA-2-GLYCOPROTEIN, IgG, IgM	BET2G	NONE	1 SST (OR PLAIN RED)	SPIN, AND SEPARATE OFF SERUM (3.0 mL)	FROZEN	NorDx - SCARBOROUGH	Immunology Lab	Tuesday
CARDIOLIPIN ANTIBODIES, IgA, IgG, IgM	CARAB	NONE						